

MEB-903

Electric Hospital Bed



Features

- Easy to clean with PP detachable platform with coating steel support (can modify to full flat surface)
- Material in reach of patient care is anti-bacterial (power coating steel)
- Full electric adjustable movement for easy handling
- Build-in control on head side rail
- ACP Nursing control panel in the end of bed
- Low voltage control box (battery for option) with battery back-up
- Removable PP head & foot board (color option: blue, wooden, dark red, yellow)
- Tuck-away half-length side rail
 - Single push button for fold down (color option: blue, wooden, dark red, yellow)
 - PP material for rail
- Gas-spring assistance
- Storage compartment
- 125 mm high quality brake caster (central lock is option)
- Urinal hooks on both side of bed (with 4 hooks)
- Auto-regression device (14 cm) of platform to prevent bed sore
- IV pole socket at 4 corners of bed
- Rotating bumper on 4 corners

Dimensions

MEB-903:

L2160 x W1010 mm (86 platform)



L2160 x W1010 mm (91 platform)

Adjustable height:

440~750 mm by PP platform

430~740 mm by steel platform

MEB-903:



Medical Master

Address: No. 59-5, Ying Cheng 1st Street, Shi-Tung Dist., Taichung City, Taiwan (ROC)

Tel: +886-4-23171028/ Fax: +886-4-23171029