

MET-6600

Transport Stretcher




Features

- Easy to clean with ABS detachable platform (with coating steel support)
- Material in reach of patient care are anti-bacterial (power coating steel)
- Full manual/ hydraulic adjustable movement for easy handling
- Gas-spring operation for back adjustment (single gas-spring)
- Hydraulic foot operation for height adjustment
- Hydraulic foot operation for Trendelenburg and reverse Trendelenburg
- Manual crank for knee adjustment
- Dual hydraulic pump (loading capacity: 350 kg)
- Foldable push handle on head & rear side
- Collapsible 3/4 length side rail
 - Single push button for fold down
 - Aluminum and stainless steel 304
- Chassis cover with store compartment and oxygen tank holder
- 200 mm high quality swivel caster
- Central locking & steering pedal on 4 sides
- 150 mm 5th wheel in the center (retractable)
- Urinal hooks on both side of bed
- Telescopic IV pole on diagonal corner (foldable)
- Rotating bumper on 4 corners
- Anti-crush bumper on both sides

Dimensions

MET-6600 Transport:
 Overall: L2100 x W800 mm
 Platform: L1870 x W620 mm
 Adjustable height: 560~910 mm
 Adjustable back: 0~88°
 Adjustable knee: 0~40°
 Trendelenburg: 0~18°
 Reverse Trendelenburg: 0~18°

MET-6600 Transport:



Medical Master

Address: No. 59-5, Ying Cheng 1st Street, Shi-Tung Dist., Taichung City, Taiwan (ROC)

Tel: +886-4-23171028/ Fax: +886-4-23171029