

# MET-6600

## Trauma Stretcher



### Features

- X-ray translucent Bakelite platform (with coating steel support) with X-ray cassette
- Full length X-ray track from head to toe (film insert from end of bed)
- Material in reach of patient care are anti-bacterial (power coating steel)
- Full manual/ hydraulic adjustable movement for easy handling
- Gas-spring operation for back adjustment (dual gas-spring)
- Hydraulic foot operation for height adjustment
- Hydraulic foot operation for Trendelenburg and reverse Trendelenburg
- Gas-spring operation for knee adjustment (only Adv. Trauma)
- Dual hydraulic pump (loading capacity: 350 kg)
- Foldable push handle on head & rear side
- Collapsible 3/4 length side rail
  - Single push button for fold down
  - Aluminum and stainless steel 304
- Chassis cover with store compartment and oxygen tank holder
- 200 mm high quality swivel caster
- Central locking & steering pedal on 4 sides
- 150 mm 5<sup>th</sup> wheel in the center (retractable)
- Urinal hooks on both side of bed
- Telescopic IV pole on diagonal corner (foldable)
- Rotating bumper on 4 corners
- Anti-crush bumper on both sides

### Dimensions

MET-6600 Trauma:

Overall: L2100 x W800 mm

Platform: L1870 x W620 mm

Adjustable height: 580~930 mm



Adjustable back: 0~88°

Adjustable knee: 0~35° (only Adv.)

Trendelenburg: 0~18°

Reverse Trendelenburg: 0~18°

MET-6600 Trauma/ Adv Trauma:



# Medical Master

Address: No. 59-5, Ying Cheng 1<sup>st</sup> Street, Shi-Tung Dist., Taichung City, Taiwan (ROC)

Tel: +886-4-23171028/ Fax: +886-4-23171029