

MEDICARE Mattress

YL-105 Mattress



- Standard size: 192 x 85/90 cm
- Thickness: 5~20 cm
- Customization size available

Polyurethane foam mattress with 30 HR density, formed by thousands of cells communicated to each other. Provides the firmness and adaptability necessary to enjoy a healthy rest, allowing a high recovery capacity.

It does not present medical contraindications, no harmful or toxic elements for health. Latex free.

It can be perfectly adjusted to articulated, hospital and geriatric beds; it does not affect in the use of the siderails. It can be provided for bariatric patients.

YL-115 Mattress



- Standard size: 192 x 85/90 cm
- Thickness: 5~20 cm
- Customization size available

The upper surface of the mattress is made of visco-elastic foam with HR55 density, with thermosensitive open cell. Due to its exceptional physical properties it maintains the memory instants after the pressure is not applied.

It perfectly adapts to the body shape and makes a mold of it, avoiding pressure points and distributing lying area, helping a correct blood circulation.

By open cell system, it allows a free air circulation and the mattress perspiration, avoiding a high level of dampness.

The lower part is made by HR30 density polyurethane foam, allowing the perfect support for the upper visco-elastic layer, assuring an optimal firmness and flexibility for the body.

Mattress cover type:

Model	YL-104	YL-115	YL-116
Type	Polyester	Polyester with PU coating	Polyester with TPU coating
Origin	Taiwan	Taiwan	Taiwan
Connection	Sewing and 1 size zipper	Welding and 3 side zipper/ sewing and 3 side zippers	Welding and 3 side zipper/ sewing and 3 side zippers
Breathable	V	V	V
Waterproof	Water Repellent	Waterproof JIS L1029:2009 AATCC 22:2100	Waterproof JIS L1029:2009 AATCC 22:2100
Fire retardant	X	CA117-2013 NFPA 701	CA117-2013 NFPA 701
Anti-bacterial	X	JIS L1902: 2015 for Staphylococcus aureus/ Klebsiella pneumonia/ E. coli	JIS L1902: 2015 for Staphylococcus aureus/ Klebsiella pneumonia/ E. coli
Elasticity	X	One side only	Two sides
Anti-slip	X	V	V
Pull belt	Option	Option	Option

Mattress foam type:

Type	MC-02	MC-04
Foam	PU foam	Top: Viscoelastic foam Bottom: PU foam
Density	30 +/-5 kg/m3	Top: 55 +/- kg/ m3 Bottom: 30 +/-5 kg/m3
Hardness	175 +/- 5 F	Top: 45 +/- F Bottom: 175 +/- 5 F
Fire retardant	X	Option
Anti-bacterial	X	Option
Elasticity	Normal	Good
Special cutting	Available	X
Egg shape	Available	X
Thickness	5~20 cm	5~20 cm